


# SERVICE

3.600 Stock version

 **Allowed user weight**  
120 kg



## *For universal application*

---

### **EFFICIENCY**

- Quick and dependable spare parts provision, secured in the long term through an expansive specialist dealer network
- Extraordinary durability, little need for repairs

### **THERAPEUTIC BENEFIT**

- For non permanent application (e.g. after an operation) as transport- and mobility aid
- Requirement based configuration through different arm- and leg support versions

### **FIELD OF USE**

- As transport- or pushing aid in hospitals, clinics and in home care
- Transport of guests with reduced mobility, e.g. in the hotel
- For short term use at home or after injuries due to an accident
- As walking aid (scutteling while sitting in a wheelchair)

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please fax the completed order form to:

05733 922 9311

Customer service 05733 922 311

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

**MEYRA**



### Basic equipment

- Drive wheel, puncture-proof, removable, quick-release axle, pressure brake
- Aluminium handrim with 15 mm distance
- Steering wheel puncture safe 200 x 50 mm
- Arm support slanted in the front with padded arm support
- Leg supports, removable and swivelling
- Incl. pedal tube
- Back straps height 400 mm
- Seat height 520 mm
- Seat depth 440 mm
- Colour: Chrome or Silverline
- Chrome version with synthetic leather covers (only i.c. with SW 450 mm)
- including calf strap, continuous



### Allowed user weight

120 kg



### Hints

Our Terms and Conditions as stated at [www.meyra.de/AGB](http://www.meyra.de/AGB) shall apply  
You can find information on the measurement process at [www.meyra.de/messverfahren](http://www.meyra.de/messverfahren)

Valid from 01.01.2024 – 31.12.2024.  
Errors, misprints and changes reserved.

When transported in public transport (e.g. train), the dimensions of the wheelchair may be higher than the guide values of the transport company. Please inquire in advance.  
When using auxiliary drives, the permissible total weight must not be exceeded.

Use our K.I.S.S.S. for orders and configurations at [www.kisss-by-meyra.de](http://www.kisss-by-meyra.de)  
You avoid queries and incorrect orders thanks to the buildability check

**Technical specifications**

Seat widths in mm	400 / 450 / 480 / 510
Width ready to drive in mm from	620 / 670 / 700 / 730
Width folded in mm from	310
Seat depth in mm	440
Seat height front/rear in mm	520 / 490
Overall height in mm	920
Back support height in mm	400
Arm support height in mm	240
Length with footrest in mm	1080
Length without footrests in mm	780
User weight	120 kg
Empty weight from	19 kg
Permissible overall weight	140 kg
Rear wheel	610 mm(24 inch)
Transport weight	11 kg
Brake test according to STANDARD	7°

Please enter your quantities in the boxes

### 01 Clinic version

Item	DESCRIPTION
9769230	Chrome/ Imitation leather in SW 450 mm

### 02 Pressure Brake

Item	DESCRIPTION
9769075	Pressure brake Silverline/Textile black in SW 400 mm
9769073	Pressure brake Silverline/Textile black in SW 450 mm
9769076	Pressure brake Silverline/Textile black in SW 480 mm
9761783	Pressure brake Silverline/Textile black in SW 510 mm

### 03 Options (delivered separately)

Item	DESCRIPTION
1072681	Height adjustable push handles Code 502 (refit kit)

### USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

### ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)

<https://www.meyra.com/gtc/>

Date

City

Signature